

# Synchronous Interferometric Measurement of Dispersion Applied to Customization of Supercontinuum Sources

H. Muñoz-Marco, J. Abreu-Afonso, V. Otgon,  
R. Dauliat, R. Jamier, P. Roy, **P. Pérez-Millán**

COST Action MP1401 WG2 Technical Meeting  
8 - 9 February, Wien, Austria

# Overview

- FYLA
- Method. Synchronous Interferometric Measurement of Dispersion
- Application. Customization of Supercontinuum Sources
- Conclusions



## FYLA

### All made in Spain, to the world

As a company, FYLA is the only one in the south of Europe, and one of the few in the world, able to produce and control the supply and manufacturing of the 100% of the components and technologies required for the production of our All-fiber pulsed Lasers, from our proprietary White Lasers Series , to pulsed lasers in the range of pico and femto seconds, with a ground breaking In-Fiber laser generation technology."

[SEE MORE](#)

### FYLA CENTRAL SITE

Ronda Guglielmo Marconi 12  
Parque Tecnológico 46980  
Paterna -Valencia -Spain  
  
Tel +34 96 136 91 90  
[fyla@fyla.com](mailto:fyla@fyla.com)

### FYLA FACTORY

Calle Formentera 24 , 08192  
Sant Quirze del Valles – Spain  
  
[fyla@fyla.com](mailto:fyla@fyla.com)



## THE COMPANY

## Motivation

- Lack of Information on D values from optical fiber manufacturers
- Limitation for labs and companies:
  - accurate modelling (pulse propagation)
  - interpretation of experimental results
  - ...



### 20/400 Ytterbium-Doped LMA Double Clad Fiber

Nufem's Large Mode Area (LMA) Ytterbium-doped double clad fiber is ideal for high power fiber lasers and amplifiers used in military, industrial, and medical applications. These fibers enable efficient, compact, diode pumped lasers that directly compete with traditional solid-state lasers. The fiber features a 20 micron diameter core and 400 micron diameter clad size with a low NA (0.08) core. With > 75% slope efficiency and compatibility with operating at > 1kW of CW output power, this fiber is ideal for use in high power single-mode industrial fiber lasers.

#### Typical Applications

- High power fiber lasers
- CW and pulsed amplifiers
- Military, industrial and medical

#### Features & Benefits

- NuCOAT™ Fluoroacrylate coating — Greater fiber durability in extreme environmental operating & storage conditions
- LMA core design — Useful for generating high CW powers
- "Fem" mode core design — Easy to maintain single mode LP01 beam through fiber & components
- PANDA-style stress structure for increased birefringence — Superior optical performance and uniformity
- All fiber proof tested to > 100 kpsi — Critical for ensuring long term reliability when coiling

#### Optical Specifications

Operating Wavelength  
Core NA  
First Cladding NA (%)  
Cladding Absorption  
Birefringence

#### PLMA-YDF-20/400-VIII

1080 – 1115 nm  
0.080 ± 0.010  
≥ 0.48  
0.50 ± 0.05 dB/m at 915 nm  
1.50 dB/m near 975 nm  
nominal  $3.5 \times 10^{-4}$

#### LMA-YDF-20/400-VIII

1080 – 1115 nm  
0.080 ± 0.010  
≥ 0.48  
0.40 ± 0.05 dB/m at 915 nm  
1.20 dB/m near 975 nm  
N/A

#### Geometrical & Mechanical Specifications

Cladding Diameter  
Cladding Diameter (Flat-to-Flat)  
Core Diameter  
Coating Diameter  
Coating Material  
Proof Test Level

400.0 ± 15.0 µm  
N/A  
20.0 ± 2.0 µm  
550.0 ± 20.0 µm  
Low Index Polymer  
≥ 100 kpsi (0.7 GN/m<sup>2</sup>)

N/A  
400.0 ± 15.0 µm  
20.0 ± 1.5 µm  
550.0 ± 20.0 µm  
Low Index Polymer  
≥ 100 kpsi (0.7 GN/m<sup>2</sup>)



A Furukawa Company

Search Website...

ESAccess™ | Contact Us | Private Site | Check Order Status | Change Region | Select another idioma ▼

[Optical Properties](#) | [Dimensions/Geometric Properties](#) | [Mechanical and Testing Data](#) | [How To Order](#) | [Typical Applications](#) | [Options](#)

#### Optical Properties

Type	Polarization-Maintaining
Cutoff Wavelength	≤ 1280 nm
Numerical Aperture (Nominal)	0.13
Mode Field Diameter	@ 1310 nm, 9.3 ± 1.0 µm
Attenuation	@ 1310 nm, ≤ 1.2 dB/km
Beat Length	@ 1310 nm, ≤ 3.8 mm
Crosstalk	@ 1310 nm/100 m, ≤ -30 dB
Crosstalk (Typical)	@ 1310 nm/100 m, ≤ -35 dB
Operating Wavelength	1310 nm

#### Dimensions/Geometric Properties

Clad Diameter	125 ± 1.0 µm
Coating/Buffer Diameter	250 ± 10 µm
Coating Type	Dual UV Acrylate
Core/Clad Offset	≤ 0.7 µm
Clad Non-Circularity	≤ 2.0%

#### Mechanical and Testing Data

Operating Temperature	-40 to +85 °C
Proof Test Level	≥ 100 kpsi 0.689 GPa



Hello guest: [Log in](#) or [Register](#) [Sales & Support](#)

[Products](#) | [Applications](#) | [Sectors](#) | [Expertise](#) | [News & Events](#) | [Corporate](#) | [Faq's](#)

#### SM FIBER

#### 125µm diameter SM specialty fibers

	SM450	SM600	SM750	SM800(5.6/125)	SM980(3.7/125)	SM980(4.5/125)	SM980(5.8/125)
Operating Wavelength (nm)	488 - 633	633 - 780	780 - 830	830 - 980	980		980 - 1500
Cut-Off Wavelength (nm)	350 - 450	500 - 600	610 - 750	660 - 800		870 - 970	
Numerical Aperture			0.10 - 0.14		0.21 - 0.23	0.17 - 0.19	0.13 - 0.15
Mode Field Diameter (µm)	2.8 - 4.1 @488nm	3.0 - 6.3 @633nm	4.5 - 6.5 @780nm	4.7 - 6.9 @830nm	3.4 - 4.0 @980nm	4.2 - 4.9 @980nm	5.3 - 6.4 @980nm
Attenuation (dB/km)	≤5.0 @488nm	≤15 @633nm	≤5.0 @780nm	≤5.0 @830nm		≤2.0 @980nm	
Proof Test (%)	1, 2 or 3 (100, 200 or 300 kpsi)						
Cladding Diameter (µm)	125 ± 1						
Core Cladding Concentricity (µm)	±0.75			±1.0			±0.50
Coating Diameter (µm)					245 ± 7		
Coating Type	Dual Acrylate						
Operating Temperature (°C)	-55 to +85						


### Current Methods

- Phase Differential Measurement of RF signals conveyed by coherent lasers
  - Spectral Resolution very high. Spectral Bandwidth very limited (eg by EOM bandwidth)
- Optical Interferometric Measurements
  - Spectral Bandwidth is (potentially) very high. Noisy, long and tedious measurements. Spectral Resolution is Low.

### Ultimate Goal

- Commercial System of Dispersion Measurement
- Full VIS-IR Band
- Universal (any optical fiber/photonic device)
- Automated
- Very Fast (minutes range)

## Method: Full VIS-IR Interferometric Measurement of Dispersion



$$I = (E_1 + E_2)(E_1 + E_2)^* = I_1 + I_2 + 2\sqrt{I_1 I_2} \cdot \cos(\phi_2(\omega) - \phi_1(\omega))$$

$$I_\omega \propto \cos(\Delta\phi(\omega))$$


$$\frac{\partial}{\partial \omega} \Delta\phi(\omega) = 0$$

$$\tau_c \cdot (L_2 - L_1) - \frac{d}{c} = \tau \cdot L$$


$$\left. \frac{\partial \beta_c}{\partial \omega} \right|_{\lambda_0} (L_2 - L_1) - \frac{d}{c} = \left. \frac{\partial \beta}{\partial \omega} \right|_{\lambda_0} L$$



$$D_c \cdot (L_2 - L_1) - \frac{1}{c} \frac{\partial d_0}{\partial \lambda} = D \cdot L$$



## Method: Full VIS-IR Interferometric Measurement of Dispersion





$$D_c \cdot (L_2 - L_1) - \frac{1}{c} \frac{\partial d_0}{\partial \lambda} = D \cdot L$$

- Experimental Curve  
 $d = d(\lambda_0)$





- Dispersion Curve



## Drawbacks

- WHITE LAMP: ALIGNMENT
- LEDs: MULTIPLE SOURCES.
- LOW SPECTRAL POWER DENSITY
- LOW VISIBILITY OF FRINGES
- NOISE
- LOW SPECTRAL RESOLUTION
- TEDIOUS
- LONG TIME CONSUMING





## Method: Full VIS-IR Interferometric Measurement of Dispersion



- Experimental Curve  
 $d = d(\lambda_0)$







- Dispersion Curve



## Advantages


- SINGLE FIBER COUPLED SOURCE
- HIGH SPECTRAL POWER DENSITY
- SYNCHRONIZED INTERFERENCES
- HIGH VISIBILITY OF FRINGES
- HIGH RESOLUTION
- LOW TIME CONSUMING
- LIABLE OF AUTOMATION

### Method: Full VIS-IR Interferometric Measurement of Dispersion



### Advantages



- SINGLE FIBER COUPLED SOURCE
- SYNCHRONIZED INTERFERENCES
- HIGH SPECTRAL POWER DENSITY
- HIGH VISIBILITY OF FRINGES
- HIGH RESOLUTION
- LOW TIME CONSUMING
- LIABLE OF AUTOMATION



# Overview

- FYLA
- Method. Synchronous Interferometric Measurement of Dispersion
- Application. Customization of Supercontinuum Sources
- Conclusions


## Fine Tuned Manufacturing of PCF to Customized SC Spectrum for OCT application





### Requirements

- Low cutoff wavelength (3dB) =  $560 \pm 10$  nm
- High cutoff wavelength (3dB) =  $850 \pm 50$  nm
- Spectral Power Density >  $200 \mu\text{W}/\text{nm}$  in range [560,850] nm

### Goal Dispersion Curve (ZDW = $967 \pm 10$ nm)




## 1. First approximation. 4 rings PCF



## 2. Fine Tuning. 7 rings PCF

Name	Material	Pressure (kPa)	Core Diameter ( $\mu\text{m}$ )	Hole Diameter ( $\mu\text{m}$ )
XLIM5	F300 silica	10.0	3.87	0.94



### Requirements **OK**

- Low cutoff wavelength (3dB) =  $560 \pm 10$  nm
- High cutoff wavelength (3dB) =  $850 \pm 50$  nm
- SPD >  $200 \mu\text{W}/\text{nm}$  in range [560,850] nm

## Outcome XLIM-FYLA Collaboration (Boosted by COST Action)

- Collaboration contract FYLA-XLIM 2017-2018
- Granted EU EUREKA Project 2018-2010
- Accepted Contribution to SPIE Photonics Europe 2018
- Several running joint experiments

# Conclusions

- Synchronous Interferometric Method to Measure Chromatic Dispersion in Optical Fibers
- Enhanced performance vs conventional methods: spectral resolution, bandwidth, simplicity, time of measurement.
- Towards universal commercial product
- Successful application to PCF manufacturing for customized SC source.

# Questions